

DharmaFECT™ siRNA Transfection Reagents

Product Description:

- Panel of four transfection reagents optimized for efficient delivery of siRNA into mammalian cells.
- Provided as a liquid at a concentration of 1 mg/mL.
- Functionally verified to produce at least 75% transfection efficiency by assessing mRNA knockdown of positive control siRNAs.
- Confirmed to be free of microbial contaminants and ribonucleases.

Shipping and Storage:

- DharmaFECT siRNA Transfection Reagents are shipped on ice packs. If the ice melts during transit there is no risk to the product. Stability testing demonstrates no loss in product functionality after prolonged storage at warm temperatures (see the Stability Testing Product Bulletin for additional details).
- Upon receipt, store at 4 °C. DO NOT FREEZE. DharmaFECT Transfection Reagents are labeled with a date of manufacture and are good for 24 months following that date when stored as recommended.

Related Products:

We recommend including a positive and negative control, such as Dharmacon RNAi Control reagents, in every RNAi experiment. For more information, click [here](#).

Accompanying Documents:

- [Basic siRNA Transfection Protocol](#).

Publication Reference Guide:

When referencing the use of DharmaFECT siRNA Transfection Reagents, please include the following information: DharmaFECT (1, 2, 3, or 4) siRNA Transfection Reagent, GE Healthcare Dharmacon, Inc., Lafayette, CO.

Product	Size	Cat. #
DharmaFECT 1	0.2 mL	T-2001-01
	0.75 mL	T-2001-02
	1.5 mL	T-2001-03
	1.5 mL × 5 tubes	T-2001-04
DharmaFECT 2	0.2 mL	T-2002-01
	0.75 mL	T-2002-02
	1.5 mL	T-2002-03
	1.5 mL × 5 tubes	T-2002-04
DharmaFECT 3	0.2 mL	T-2003-01
	0.75 mL	T-2003-02
	1.5 mL	T-2003-03
	1.5 mL × 5 tubes	T-2003-04
DharmaFECT 4	0.2 mL	T-2004-01
	0.75 mL	T-2004-02
	1.5 mL	T-2004-03
	1.5 mL × 5 tubes	T-2004-04
DharmaFECT 1, 2, 3, and 4	0.2 mL each	T-2005-01
	0.75 mL each	T-2005-02
	1.5 mL each	T-2005-03

GE, imagination at work and GE monogram are trademarks of General Electric Company. Dharmacon, Inc., a General Electric company doing business as GE Healthcare. All other trademarks are the property of General Electric Company or one of its subsidiaries. ©2014 General Electric Company—All rights reserved. First version published June 2014. GE Healthcare UK Limited, Amersham Place, Little Chalfont, Buckinghamshire, HP7 9NA, UK



Orders can be placed at:
gelifesciences.com/dharmacon

Customer Support: cs.dharmacon@ge.com
Technical Support: ts.dharmacon@ge.com or
1.800.235.9880; 303.604.9499 if you have any questions.